






corpuls **cpr** Daily Function Checks

Date: _____

Step	Task	Result	Check Box*
1	Visually check the corpuls cpr and its accessories for damage	<ul style="list-style-type: none"> The corpuls cpr and its accessories are not visibly damaged 	
2	Check the charge status of the battery and reserve battery	<ul style="list-style-type: none"> At least 4 LEDs of the battery display should illuminate Charge the battery and / reserve battery if indicated 	
3	Check whether an SD card is inserted	<ul style="list-style-type: none"> An SD Card is inserted Insert an SD card if the SD card housing is empty 	
4	Assemble the arm (with the battery inserted) on the board	<ul style="list-style-type: none"> The arm locks easily into the socket of the board The arm is securely connected to the board 	
5	Release the locking lever	<ul style="list-style-type: none"> The arm is easy to move and does not lock 	
6	Close the locking lever	<ul style="list-style-type: none"> All moveable parts of the arm are fixed 	
7	Check the stamp for stress whitening and insert stamp**	<ul style="list-style-type: none"> The stamp clicks into position 	
8	Power on the arm	<ul style="list-style-type: none"> The self-test starts The screen is illuminated The Start / Stop button flashes white and a short audible signal is emitted The arm does not signal any alarm The main screen or start screen appears 	
9	Perform stamp position check		
	a) Open the locking lever	<ul style="list-style-type: none"> The LED of the Start / Stop button lights up yellow The display shows the symbol  	
	b) Press the stamp gently inwards	<ul style="list-style-type: none"> The LED of the Start / Stop button lights up green The display shows no symbol 	
	c) Press the stamp firmly inwards	<ul style="list-style-type: none"> The LED of the Start / Stop button lights up red The display shows the symbol  	
	d) Release the stamp	<ul style="list-style-type: none"> The LED of the Start / Stop button lights up yellow The display shows the symbol  	
	e) Start therapy	<ul style="list-style-type: none"> The therapy does not start The display shows the symbol  	
	f) Close the locking lever	<ul style="list-style-type: none"> The symbol  disappears from the display 	
g) Start therapy	<ul style="list-style-type: none"> The arm performs compressions in the configured mode *** 		
10	Power off the arm	<ul style="list-style-type: none"> The arm does not signal any alarm The arm powers off 	

*Laminate this document to use the **Check Box** when completing the individual components of the **Daily Function Checks** for the **corpuls cpr**.

**The stamp is intended for repeated use; The user must replace the stamp immediately if stress whitening occurs. Stress whitening are fine cracks in the material which typically occur on the stamp surface. If the user does not replace the stamp in case of stress whitening, the stamp can break along these cracks and injure the patient.

***Thorax compression devices can caused various operation noises due to differences of compression rate levels. Squeaking noises are in general normal operation noises of the **corpuls cpr**.

