MEDUMAT Standard²

Ventilator

Safety information for devices with CO₂ measurement
Supplement to the instructions for use

This is a supplement to the instructions for use. Please read this safety information along with the instructions for use before using the product. Failure to observe the information contained therein can result in serious injuries or even death.

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Safety information

This safety information supplements the “Safety information” chapter of the instructions for use for the device. Please also always observe the instructions for use. The instructions for use and this supplement can be found at: www.weinmann-emergency.com.

 **WARNING****Delay in treatment due to abnormal CO₂ measurement during chest compressions!**

The device's CO₂ measurement function is **not** designed for the high frequencies (100/min to 120/min) encountered during chest compressions. Chest compression can be detected by the device as high-frequency breathing, which can result in the etCO₂ measurement being erroneously displayed as too low and trigger implausible alarms. This can confuse the user and delay treatment. ⇒ Do **not** rely on the etCO₂ measurement as the sole indication for stopping cardiopulmonary resuscitation!

Please note: The etCO₂ measurement and the CO₂ curve can still be used for:

- checking the tube position;
- identifying the return of spontaneous circulation (ROSC) by a sudden increase in the etCO₂ measurement.

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Medical Technology GmbH + Co. KG