






corpuls **cpr** Alarms Quick Guide

High Priority Alarms




Must be actioned immediately. Therapy will stop. Therapy settings cannot be adjusted during an active or suspended alarm. The alarm must be dealt with prior to restarting the therapy.

Common Alarms

Alarm Symbol	Cause	Outcome	Required Action
	The locking lever was not initially opened after the arm was switched on	<ul style="list-style-type: none"> Therapy not possible 	<ul style="list-style-type: none"> Open the locking lever Adjust the arm and set the therapy position Close the locking lever Start therapy
	Locking lever not closed	<ul style="list-style-type: none"> Therapy not possible 	<ul style="list-style-type: none"> Close the locking lever
	Arm set too low	<ul style="list-style-type: none"> Therapy not possible 	<ul style="list-style-type: none"> Raise the arm Check the stamp position and correct
	Therapy stopped for longer than 8 seconds	<ul style="list-style-type: none"> No therapy performed on the patient 	<ul style="list-style-type: none"> Continue therapy with the corpulscpr Confirm the alarm OR Continue CPR manually without the corpulscpr
	The battery charge status is critically low	<ul style="list-style-type: none"> Imminent device power off will occur 	<ul style="list-style-type: none"> Insert a fully charged reserve battery OR Apply external power source OR Continue CPR manually without the corpulscpr

High Priority Alarms




Less Common Alarms

Alarm Symbol	Cause	Outcome	Required Action
	Bluetooth connection to corpuls3 or corpuls3T has been lost	<ul style="list-style-type: none"> Synchronised mode lost 	<ul style="list-style-type: none"> Manually start / stop the therapy at the corpulscpr
	Temperature of the arm is critically high	<ul style="list-style-type: none"> Arm temperature exceeds limitations Imminent device power off will occur 	<ul style="list-style-type: none"> Continue CPR manually without the corpulscpr <p><small>⚠ During the power down procedure, the alarm symbol will re appear as a message.</small></p>
	System error	<ul style="list-style-type: none"> Injury of the patient possible Alarm failure Sensor failure Therapy failure 	<ul style="list-style-type: none"> Switch off the corpulscpr Continue CPR manually without the corpulscpr Contact customer service <p><small>⚠ During the power down procedure, the alarm symbol will re appear as a message.</small></p>

Medium Priority Alarms





Therapy may continue during a medium priority alarm, however the alarm must be actioned quickly.

Common Alarms

Alarm Symbol	Cause	Outcome	Required Action
	Arm too high	<ul style="list-style-type: none"> The configured compression depth cannot be achieved Inadequate therapy 	<ul style="list-style-type: none"> Open the locking lever Adjust the arm and reset the therapy position Close the locking lever Start therapy
	Therapy stopped	<ul style="list-style-type: none"> No therapy performed on the patient The alarm escalates to a HIGH priority alarm after 8 seconds 	<ul style="list-style-type: none"> Continue therapy with the corpulscpr Confirm the alarm OR Continue CPR manually without the corpulscpr
	Battery charge status is low	<ul style="list-style-type: none"> Device power off will occur 	<ul style="list-style-type: none"> Insert a fully charged reserve battery OR Apply external power source OR Continue CPR manually without the corpulscpr

Medium Priority Alarms

Less Common Alarms

Alarm Symbol	Cause	Outcome	Required Action
	Overload	<ul style="list-style-type: none"> The corpulscpr arm is possibly overloaded 	<ul style="list-style-type: none"> Check the position of the stamp Adjust the arm and reset the therapy position Confirm the alarm <p> Consider that the rigidity of the patient's thorax may be higher than normal.</p>
	Temperature of the arm is too high	<ul style="list-style-type: none"> Arm temperature exceeds limitations 	<ul style="list-style-type: none"> Adjust the compression depth and compression rate If necessary, continue manual CPR without the corpulscpr
	The start-up self-test has failed	<ul style="list-style-type: none"> The technical specifications of the device cannot be ensured 	<ul style="list-style-type: none"> Do not use the corpulscpr Continue manual CPR without the corpulscpr Contact customer service

 **Remove the battery if the cpr arm cannot be powered off in a critical malfunction emergency** 