


C3T Daily Function Checks


The function check comprises a visual inspection of the outer housing and a function check of the functionalities and options of the **C3T**. All the modules of the **corpuls3T** must be connected mechanically.

Additionally, inspect the test box and cable / test load for visible damage prior to use.


Do not use the testload / testbox if visibly damaged.


Powering On


 Press the **On / Off** key to power the device in compact mode


 *An Internal test will be automatically complete*

Module Connection


 Separate the modules


 *Observe that the connection status symbol changes from a bar to a wave symbol*

 Reconnect the modules

 *Observe that the connection status symbol changes from a wave to a bar symbol*

Defibrillator / Pacer


-  Connect the testload or testbox to the therapy socket
- ⇒ Press the [Yes] softkey to confirm the message **"Mark mission as test mission?"**
 - ⇒ Select Manual mode
 - Select an energy level of 50 J for the testload
 - Or 200 J for the testbox
 - ⇒ Charge
 - ⇒ Deliver the shock


 *Check the shock protocol printout for the message*


"Test successful" OR "Test failed"


The energy is internally discharged via the testload / testbox


Printer

 Ensure sufficient printer paper is available and the end of roll markings (red stripe) are not visible

 *Reload a new roll if the red stripe is visible*


 Run a test printout, e.g. a screenshot


 *The printout is legible and of good quality*


 Ensure a replacement roll is present


 *A replacement roll is present*


Power Supply


 Check state of charge of the compact device when powered by the batteries only


 *State of charge of the battery is greater than 30%*


 Connect to mains supply

 Check state of charge of the compact device, observed in %

 Separate all 3 modules while powered on & Remove the device from mains power supply





 Ensure each module has a battery inserted

 Check the remaining run time of the monitoring unit and the patient box


 *Remaining runtime of the modules > 120 minutes.*



Accessory Checks

-  ✓ All patient monitoring and therapy intermediate cables are present, plugged into the device and undamaged.
 -  ✓ There is an appropriate supply of all patient monitoring and therapy accessories.
 -  ✓ All re-usable accessories are undamaged.
 -  ✓ Packaging for single patient use consumables is intact, undamaged and within the expiry date.
-

CO₂ Function checks

-  With the **C3T** powered on, function test the CO₂ cable and sensor:
 - ⇒ Disconnect the CO₂ sensor from its transport adaptor
 - ⇒ Observe for the immediate visual and audible alert **“CO2 sensor loose”**, indicating that the CO₂ cable and sensor are in working order.
-

Powering Off

-  Press the **On / Off** key to power off the device in compact mode

If the option **“Mission upload”** is enabled and configured:

- ⇒ Press the **[OK]** softkey to start the upload of mission recorded on the CF card to the configured server
- ⇒ Allow sufficient time for the data upload to be completed
- ⇒ The **C3T** will power off on completion of the data upload

Avoid manually powering off the device while the mission data upload is taking place!

