






corpuls **cpr** Daily Function Checks

Date:

Step	Task	Result	Check Box*
1	Check the corpuls cpr and its accessories for damage	<ul style="list-style-type: none"> The corpuls cpr and its accessories are not visibly damaged 	
2	Check the charge status of the battery and spare battery	<ul style="list-style-type: none"> At least four LEDs of the battery display should illuminate Charge depleted batteries 	
3	Assemble the arm (with battery inserted) on the board	<ul style="list-style-type: none"> The arm locks easily into the socket of the board The arm is securely connected to the board 	
4	Check the stamp for stress whitening and insert stamp**	<ul style="list-style-type: none"> The stamp clicks into position 	
5	Power on the arm	Device self-test is performed: <ul style="list-style-type: none"> Display screen is illuminated The therapy button flashes white, and a short audible signal is emitted No alarm is signalled The main screen or start screen appears 	
6	Check SD card is present (device will alarm if card is missing)	<ul style="list-style-type: none"> SD card is present Insert SD card if necessary 	
7	Perform stamp position check		
	a) Open the locking lever	<ul style="list-style-type: none"> The arm moves easily up & down, horizontally, and flexes at the pivot joint The LED of the therapy button lights up yellow The display shows the symbol  	
	b) Press the stamp gently inwards	<ul style="list-style-type: none"> The LED of the therapy button lights up green The display shows no symbol 	
	c) Press the stamp firmly inwards	<ul style="list-style-type: none"> The LED of the therapy button lights up red The display shows the symbol  	
	d) Release the stamp	<ul style="list-style-type: none"> The LED of the therapy button lights up yellow The display shows the symbol  	
	e) Start therapy	<ul style="list-style-type: none"> Therapy does not start The display shows the symbol  	
	f) Close the locking lever	<ul style="list-style-type: none"> All movement is locked The symbol  disappears from the display 	
	g) Start therapy	<ul style="list-style-type: none"> The arm performs compressions in the configured mode 	
8	Power off the arm	<ul style="list-style-type: none"> The device does not signal any alarm The device powers off 	

*Laminate this document to use the **Check Box** when completing the individual components of the **Daily Function Checks** for the **corpuls cpr**.

**The stamp is intended for repeated use; The user must replace the stamp immediately if stress whitening occurs. Stress whitening are fine cracks in the material which typically occur on the stamp surface. If the user does not replace the stamp in case of stress whitening, the stamp can break along these cracks and injure the patient.



Ortus Academy

Ortus Academy is a trading name of Ortus Technology Limited.
Company Registration 07554900

ORTUS
ACADEMY

Academy corpuls cpr_Daily Function Checks_SWV1.3_A008_v1

This document has been developed to provide users with an overview for safe device operation in accordance with its function and intended use and does not replace the user manual. All persons dealing with use, maintenance and troubleshooting must read and implement the user manual. Illustrations may not reflect that displayed on your device. Subject to technical modifications, mistakes, and printing errors.