



corpuls | GS Elektromedizinische Geräte G. Stemple GmbH | Hauswiesenstr. 26 | 86916 Kaufering

TO WHOM IT MAY CONCERN

Ihr Zeichen • Your reference    Unser Zeichen • Our reference    Direktwahl • Direct dialing    E-Mail • e-mail    Datum • Date  
15. December 2023

Dear Ladies and Gentlemen,

with this letter we hereby confirm that the following list of cleaning and disinfecting agents were verified by GS Elektromedizinische Geräte G. Stemple GmbH to be compatible with the materials of the corpuls cpr devices. Therefore, they can be applied safely to the corpuls cpr for reprocessing purposes.

Other cleaning and disinfecting agents may be used. In that case, GS corpuls cannot take responsibility for possible material damage of the devices. In order to minimize the risk of such damage, we recommend to choose agents that contain the same active ingredients in a concentration not exceeding the ones from the list below and to perform your own material compatibility verification.

- 1. Ultrasol® Oxywipes
- 2. Clinell® Universalwipes
- 3. Mikrobac® Forte
- 4. Kohrsolin® FF
- 5. Dismozon® plus
- 6. Bacillol® Sensitive 30
- 7. Mikrobac® Virucidal Tissues

Best regards,

---

Dr.-Ing. Lukas von Stülpnagel  
Head of Product Management