



corpuls | GS Elektromedizinische Geräte G. Stemple GmbH | Hauswiesenstr. 26 | 86916 Kaufering

## TO WHOM IT MAY CONCERN

Ihr Zeichen • Your reference

Unser Zeichen • Our reference

Direktwahl • Direct dialing

E-Mail • e-mail

Datum • Date

15. December 2023

Dear Ladies and Gentlemen,

with this letter we hereby confirm that the following list of cleaning and disinfecting agents were verified by GS Elektromedizinische Geräte G. Stemple GmbH to be compatible with the materials of the corpuls3 devices. Therefore, they can be applied safely to the corpuls3 for reprocessing purposes.

1. Ultrasol® Oxywipes
2. Clinell® Universalwipes
3. Mikrobac® Forte
4. Kohrsolin® FF
5. Dismozon® plus
6. Bacillol® Sensitive 30
7. Mikrobac® Virucidal Tissues

Best regards,

---

Dr.-Ing. Lukas von Stülpnagel  
Head of Product Management