



CPR Feedback Change Notice




Yorkshire Ambulance Service

There will be 2 technologies of CPR feedback cables and sensors circulating in the YAS fleet of **C3**'s.
Below are the differences between the cables and CPR feedback sensors for the **current** and **new** technology setup.
They are **not interchangeable**.

Current Technology

corPatch CPR Feedback

Identified by **BLUE** connectors



New Technology

PrimeCPR Feedback

Identified by **GREEN** connectors



Ortus Academy

Ortus Academy is a trading name of Ortus Technology Limited.
Company Registration 07554900

ORTUS
ACADEMY

YAS_C3_CPR Feedback_SPU_Change Notice_SWV4.3_Z003_v1

This document has been developed to provide users with an overview for safe device operation in accordance with its function and intended use and does not replace the user manual. All persons dealing with use, maintenance and troubleshooting must read and implement the user manual. Illustrations may not reflect that displayed on your device. Subject to technical modifications, mistakes, and printing errors.