

corpuls **cpr** Arm

Standard Configuration

This document details the **Standard Configuration** settings applied to the corpuls **cpr** arm.

Settings	Standard Configuration
Mode	30:2
Compression depth	5.5 cm
Compression rate	100 1/min
Usage selection	Advanced
User level	Default
Duration of pause for ventilation	4 s
Volume	10
LED brightness	7
Backlight	7
Start screens	Not Active
Key tone	Enabled
Ventilation signal frequency	10 1/ min
Duration of alarm suspension	120 s
Code for user level DEFAULT	0000
Code for user level OPERATOR	1234
Bluetooth	Active
BT PIN	1234
LED	Enabled
Audio (Tone)	Enabled
Reminder signal	Enabled
Overload alarm	Disabled
Language	English
Date	YYYY-MM-DD
Time	HH:MM:SS